

Sample Name :Melittin,honey bee

Cat No. : RP20415

Lot No. : E8006423D

Time Processed :11:08:34

Month-Day-Year Processed :04/25/2023

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	Pump A B.Conc	5
10.00	Pumps	Pump A B.Conc	35
35.00	Pumps	Pump A B.Conc	95
37.00	Pumps	Pump A B.Conc	95
37.01	Pumps	Pump A B.Conc	5
45.00	Pumps	Pump A B.Conc	5
45.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: SS-CM-0286

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	21.796	14902	5100	0.136
2	21.960	10864743	1254614	98.831
3	22.279	82163	8608	0.747
4	23.188	31420	2571	0.286
Total		10993229	1270893	100.000