

Sample Name:TAT peptide
Cat No. : RP20256
Lot No. : E8108226G
Time Processed :14:25:12
Month-Day-Year Processed :07/11/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm
<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm
Equipment: ZJ19010141

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	4.071	47900	5134	0.381
2	4.299	210258	19131	1.674
3	4.729	238186	43309	1.896
4	5.218	11945298	744918	95.098
5	5.456	49318	9605	0.393
6	5.647	20884	2965	0.166
7	6.058	49258	7931	0.392
Total		12561102	832991	100.000