


Mass Spectrum



Sample Information

Month-Day Processed : 07/11/26
Time Processed : 14:02:49
Injection Volume : 0.3
Sample Name : TAT peptide

Cat No. : RP20256

Lot No. : E8108226G

Theoretical MW : 1621.91

Observed MW : 1621.6

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : ZJ21010035
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH