

Sample Name :OVA Peptide (323-339)
Cat No. :RP10610
Lot No. :E8170325L
Time Processed :14:52:31
Month-Day-Year Processed :12/27/2025
Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm
<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ25010053

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	8.308	6923114	1055974	95.110
2	8.575	217509	25395	2.988
3	8.725	100296	13656	1.378
4	8.950	15663	2365	0.215
5	9.017	6881	1824	0.095
6	10.600	15579	1929	0.214
Total		7279043	1101143	100.000