


Mass Spectrum



Sample Information

Month-Day Processed : 01/18/23
Time Processed : 19:50:30
Injection Volume : 0.2
Sample Name : TAT peptide
Cat No. : RP20256
Lot No. : E8259123A
Theoretical MW : 1621.91
Observed MW : 1622.0

Interface :ESI
Nebulizing Gas Flow :1.5L/min
CDL Temp :250
Block Temp :200

Equipment :ZJ22010150
Interface Bias :+4.5 kV
Drying Gas Flow :5 L/min
T.Flow :0.2 ml/min
B.conc :50%H2O/50%MeOH