

Sample Name :Influenza A NP(366-374) Strain A/PR/8/35

Cat No. :RP20267

Lot No. :E8270226G

Time Processed :20:04:45

Month-Dav-Year Processed :07/20/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: ZJ21010375

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	8.758	30615	5757	1.431
2	9.017	3720	459	0.174
3	9.208	4683	651	0.219
4	9.375	4057	671	0.190
5	9.683	2056846	367476	96.133
6	9.933	39665	4531	1.854
Total		2139586	379544	100.000