

Sample Name :Melittin,honey bee
Cat No. : RP20415
Lot No. : E8305323G
Time Processed :4:30:29
Month-Day-Year Processed :07/06/2023
Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: GK12010012

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	9.357	12794	537	0.206
2	18.884	84439	7631	1.362
3	19.691	49383	5618	0.796
4	20.453	6038168	575906	97.368
5	21.786	16586	1793	0.267
Total		6201370	591485	100.000