


Mass Spectrum



Sample Information		Equipment	
Month-Day Processed :	07/06/23	Interface	: ZJ21010035
Time Processed :	9:00:36	Interface Bias	: +4.5 kV
Injection Volume :	0.3	Drying Gas Flow	: 5 L/min
Sample Name :	Melittin,honey bee	T.Flow	: 0.2 ml/min
Cat No. :	RP20415	B.conc	: 50%H2O/50%MeOH
Lot No. :	E8305323G		
Theoretical MW :	2847.50		
Observed MW :	2847.6		