

Sample Name :Glucagon-Like Peptide (GLP) II, human

Cat No. : RP10769

Lot No. : E8311724G

Time Processed :0:38:48

Month-Day-Year Processed :07/24/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: SS-CM-0309

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	17.850	61759	5645	0.290
2	19.417	67666	7414	0.317
3	19.925	188275	14926	0.883
4	20.317	382511	37334	1.794
5	20.625	154915	18724	0.727
6	20.925	20340291	1803645	95.421
7	21.433	48433	5690	0.227
8	21.658	72405	5279	0.340
Total		21316256	1898656	100.000