


Mass Spectrum



Sample Information

|                       |                                       |                     |           |                 |                 |
|-----------------------|---------------------------------------|---------------------|-----------|-----------------|-----------------|
| Month-Day Processed : | 07/24/24                              | Interface           | :ESI      | Equipment       | : ZJ21010035    |
| Time Processed :      | 10:35:24                              | Nebulizing Gas Flow | :1.5L/min | Interface Bias  | : +4.5 kV       |
| Injection Volume :    | 0.4                                   | CDL Temp            | :250      | Drying Gas Flow | :5 L/min        |
| Sample Name :         | Glucagon-Like Peptide (GLP) II, human | Block Temp          | :200      | T.Flow          | :0.2 ml/min     |
| Cat No. :             | RP10769                               |                     |           | B.conc          | :50%H2O/50%MeOH |
| Lot No. :             | E8311724G                             |                     |           |                 |                 |
| Theoretical MW :      | 3766.20                               |                     |           |                 |                 |
| Observed MW :         | 3765.9                                |                     |           |                 |                 |