

Sample Name :OVA Peptide (323-339)

Cat No. :RP10610

Lot No. :E8318426H

Time Processed :13:18:45

Month-Day-Year Processed :08/10/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	Solvent B Conc.	5
25.00	Pumps	Solvent B Conc.	65
25.01	Pumps	Solvent B Conc.	95
27.00	Pumps	Solvent B Conc.	95
27.01	Pumps	Solvent B Conc.	5
33.00	Pumps	Solvent B Conc.	5
33.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ19010015

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	7.242	19602944	1640369	97.353
2	7.692	253194	32111	1.257
3	7.908	60019	8452	0.298
4	8.033	83491	7191	0.415
5	10.983	45917	4821	0.228
6	11.308	54935	8376	0.273
7	12.867	35530	5133	0.176
Total		20136031	1706453	100.000