

Sample Name :Glucagon-Like Peptide (GLP) II, human
Cat No. : RP10769
Lot No. : E8319526H

Time Processed :12:49:23
Month-Day-Year Processed :08/09/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm
Equipment: ZJ25010054

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	18.100	13452	2220	0.243
2	18.317	11234	2024	0.203
3	18.558	65144	8382	1.177
4	18.908	36924	5128	0.667
5	19.092	35873	8769	0.648
6	19.283	5272433	613772	95.247
7	19.767	18234	2421	0.329
8	19.975	82226	7878	1.485
Total		5535518	650594	100.000