

Sample Name :OVA Peptide (323-339)

Cat No. : RP10610

Lot No. : E8391226E

Time Processed :10:44:35 AM

Month-Day-Year Processed :05/17/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ17010508

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	9.900	15467	2493	0.212
2	10.058	98375	13494	1.347
3	10.375	6955193	916810	95.235
4	10.575	176615	42712	2.418
5	10.733	16701	3914	0.229
6	10.817	12420	2736	0.170
7	11.058	28386	3022	0.389
Total		7303158	985181	100.000