

Sample Name:KRAS G12V Peptide (VVVGAVGVGK)

Cat No. : RP30229

Lot No. : E8403226E

Time Processed :20:30:07

Month-Day-Year Processed :05/18/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ19010324

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	7.125	59312	8406	0.464
2	7.325	44129	5242	0.346
3	7.608	18020	3461	0.141
4	7.842	12266143	796262	96.062
5	8.375	96246	10587	0.754
6	9.108	9392	1756	0.074
7	9.167	15682	2032	0.123
8	9.450	260121	29969	2.037
Total		12769045	857715	100.000