


Mass Spectrum



Sample Information

Month-Day Processed :	03/25/23	Equipment	:ZJ22010150
Time Processed :	8:59:36	Interface Bias	: +4.5 kV
Injection Volume :	0.1	Drying Gas Flow	:5 L/min
Sample Name :	Amylin, human, amide	T.Flow	:0.2 ml/min
Cat No. :	RP11278	B.conc	:50%H2O/50%MeOH
Lot No. :	E8403523C		
Theoretical MW:	3903.28		
Observed MW :	3903.9		