


Mass Spectrum



Sample Information

Month-Day Processed : 04/04/24
Time Processed : 10:33:49
Injection Volume : 0.1
Sample Name : Exendin-4
Cat No. : RP10874
Lot No. : E8440024C
Theoretical MW: 4186.66
Observed MW : 4186.0

Interface

Nebulizing Gas Flow :1.5L/min
CDL Temp :250
Block Temp :200

Equipment

Interface Bias
Drying Gas Flow
T.Flow
B.conc

:ZJ22010150

: +4.5 kV
:5 L/min
:0.2 ml/min
:50%H2O/50%MeOH