


Mass Spectrum



Sample Information

Month-Day Processed : 03/19/23
Time Processed : 18:37:31
Injection Volume : 0.2
Sample Name : TAT peptide
Cat No. : RP20256
Lot No. : E8464823C

Theoretical MW:1621.91
Observed MW :1622.0

Interface :ESI
Nebulizing Gas Flow :1.5L/min
CDL Temp :250
Block Temp :200
Equipment : ZJ21010035
Interface Bias : +4.5 kV
Drying Gas Flow :5 L/min
T.Flow :0.2 ml/min
B.conc :50%H₂O/50%MeOH