

Sample Name :Calcitonin Gene Related Peptide (CGRP) (8-37), rat
Cat No. : RP11090
Lot No. : E8581526G
Time Processed :6:44:50
Month-Day-Year Processed :07/04/2026


Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>
<Detector A>
Column :Inertsil ODS-3 4.6 x 250 mm
Equipment: ZJ19010325

<Chromatogram>



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	10.809	13310476	1282532	98.886
2	11.104	98240	11300	0.730
3	11.402	24199	2313	0.180
4	12.350	27532	930	0.205
Total		13460446	1297076	100.000