

Sample Name :C-Type Natriuretic Peptide (CNP) (1-22), human

Cat No. : RP11110

Lot No. : E8761324J

Time Processed :6:24:31

Month-Day-Year Processed :11/05/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: SS-CM-0309

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	14.027	60891	10195	0.489
2	14.773	12064108	1368134	96.978
3	15.053	27205	5174	0.219
4	15.627	27909	3807	0.224
5	16.747	259990	30814	2.090
Total		12440104	1418125	100.000