


Mass Spectrum



Sample Information		Equipment : ZJ21010035	
Month-Day Processed :	11/05/24	Interface Bias	: +4.5 kV
Time Processed :	12:14:57	Drying Gas Flow	: 5 L/min
Injection Volume :	0.3	T.Flow	: 0.2 ml/min
Sample Name :	C-Type Natriuretic Peptide (CNP) (1-22), human	B.conc	: 50%H2O/50%MeOH
Cat No. :	RP11110		
Lot No. :	E8761324J		
Theoretical MW:	2197.60		
Observed MW :	2197.5		