

Sample Name : Glucagon-Like Peptide (GLP) II, human
 Cat No. : RP10769
 Lot No. : E8844322A

Time Processed : 19:29:56
 Month-Day-Year Processed : 01/29/2022


Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow: 1 ml/min
 Wavelength: 220 nm

Time	Module	Command	Value
0.01	Pumps	Pump A B.Conc	5
25.00	Pumps	Pump A B.Conc	65
25.01	Pumps	Pump A B.Conc	95
27.00	Pumps	Pump A B.Conc	95
27.01	Pumps	Pump A B.Conc	5
35.00	Pumps	Pump A B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm



1 Detector A Channel 1 / 220nm

Peak Table

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	13.237	1853	297	0.081
2	13.371	4233	638	0.184
3	13.550	9084	1371	0.395
4	14.196	34015	3416	1.480
5	14.395	2240895	301394	97.496
6	15.393	1869	294	0.081
7	15.996	2655	286	0.116
8	16.208	3854	485	0.168
Total		2298457	308182	100.000