


Mass Spectrum



Sample Information

Month-Day Processed :	01/29/22	Equipment	: GK11010007
Time Processed :	10:08:47	Interface Bias	: +4.5 kV
Injection Volume :	0.2	Drying Gas Flow	: 5 L/min
Sample Name :	Glucagon-Like Peptide (GLP) II, human	T.Flow	: 0.2 ml/min
SCat No. :	RP10769	B.conc	: 50% H ₂ O/50% MeOH
Lot No. :	E8844322A		
Theoretical MW :	3766.20		
Observed MW :	3766.0		