

Sample Name :TAT peptide
 Cat No. : RP20256
 Lot No. : E8863323L
 Time Processed :2:49:39
 Month-Day-Year Processed :12/20/2023
 Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm
 Equipment: ZJ20010139

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	5.800	203557	28362	1.846
2	6.142	10708615	915935	97.130
3	12.642	70343	8793	0.638
4	15.033	6234	1057	0.057
5	15.417	9523	1250	0.086
6	15.558	7400	1165	0.067
7	15.692	19308	1785	0.175
Total		11024978	958346	100.000