


Mass Spectrum



Sample Information		Equipment	: ZJ21010035
Month-Day Processed :	12/19/23	Interface Bias	: +4.5 kV
Time Processed :	17:06:34	Drying Gas Flow	: 5 L/min
Injection Volume :	0.3	T.Flow	: 0.2 ml/min
Sample Name :	TAT peptide	B.conc	: 50% H ₂ O/50% MeOH
Cat No. :	RP20256		
Lot No. :	E8863323L		
Theoretical MW :	1621.91		
Observed MW :	1622.0		