

Sample Name :Melittin,honey bee
Cat No. : RP20415
Lot No. : E8863725E
Time Processed :7:27:11
Month-Day-Year Processed :05/25/2025
Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ21010376

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	16.583	19526	1595	0.206
2	18.325	91415	5909	0.967
3	19.017	51694	2600	0.547
4	19.333	167318	9225	1.769
5	19.925	9127053	727485	96.511
Total		9457006	746813	100.000