

Sample Name :Influenza A NP(366-374) Strain A/PR/8/35
 Cat No. : RP20267
 Lot No. : E8988726F
 Time Processed :22:39:17
 Month-Day-Year Processed :06/18/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	Solvent B Conc.	5
25.00	Pumps	Solvent B Conc.	65
25.01	Pumps	Solvent B Conc.	95
27.00	Pumps	Solvent B Conc.	95
27.01	Pumps	Solvent B Conc.	5
33.00	Pumps	Solvent B Conc.	5
33.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm
 Equipment: ZJ19010015

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	6.017	2462	898	0.045
2	6.242	5467873	664465	99.163
3	6.508	29831	5844	0.541
4	6.908	7315	438	0.133
5	7.417	6537	1062	0.119
Total		5514017	672706	100.000