

Sample Name :Biotin-β-Amyloid (1-40)

Cat No. : RP30249

Lot No. : E9007026D

Time Processed :21:01:54

Month-Day-Year Processed :04/17/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :UniPSN 30-300 30um, 4.6x250mm

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	14.483	231099	6100	0.379
2	16.058	2176134	77173	3.568
3	16.992	57958690	1206407	95.041
4	18.658	617146	17320	1.012
Total		60983068	1307001	100.000