


Mass Spectrum



Sample Information
Month-Day Processed : 04/10/22
Time Processed : 7:32:59
Injection Volume : 0.2
Sample Name : C-Type Natriuretic Peptide (CNP) (1-22), human
Cat No. : RP11110
Lot No. : E9029622D
Theoretical MW : 2197.6
Observed MW : 2197.5

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : ZJ21010035
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%M