


Mass Spectrum



Sample Information

Month-Day Processed :	03/27/23	Equipment	: ZJ21010035
Time Processed :	17:11:16	Interface Bias	: +4.5 kV
Injection Volume :	0.4	Drying Gas Flow	: 5 L/min
Sample Name :	KRAS G12V Peptide (VVVGAVGVGK)	T.Flow	: 0.2 ml/min
Cat No. :	RP30229	B.conc	: 50% H ₂ O/50% MeOH
Lot No. :	E9073023C		

Theoretical MW: 884.08

Observed MW : 884.0