


Mass Spectrum



Sample Information

Month-Day Processed : 03/27/23  
Time Processed : 17:56:08  
Injection Volume : 0.4  
Sample Name : SARS-CoV-2 Spike-derived Peptide (NGVEGFNCYFPPLQSY)  
Cat No. : RP30234  
Lot No. : E9123523C  
Theoretical MW : 1737.89  
Observed MW : 1737.6

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200  
Equipment : ZJ21010035  
Interface Bias : -3.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50%H2O/50%MeOH