


Mass Spectrum



Sample Information

Month-Day Processed :	03/27/23	Equipment	:ZJ22010150
Time Processed :	17:45:37	Interface Bias	: +4.5 kV
Injection Volume :	0.1	Drying Gas Flow	:5 L/min
Sample Name :	KRAS WT Peptide (VVVGAGGVGK)	T.Flow	:0.2 ml/min
Cat No. :	RP30235	B.conc	:50%H2O/50%MeOH
Lot No. :	E9126823C		
Theoretical MW:842.00			
Observed MW : 842.2			