

Sample Name:Glucagon-Like Peptide (GLP) II, human

Cat No. : RP10769

Lot No. : E9174425G

Time Processed :20:20:17

Month-Day-Year Processed :07/10/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ20010140

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	14.692	9700	1663	0.151
2	18.050	35212	5109	0.549
3	18.550	10814	2857	0.169
4	18.633	38662	5145	0.603
5	18.917	129728	16439	2.024
6	19.342	19563	3972	0.305
7	19.683	6098344	756365	95.156
8	20.183	35627	4401	0.556
9	20.425	21865	2747	0.341
10	20.650	9264	1460	0.145
Total		6408778	800158	100.000