

Sample Name:C-Type Natriuretic Peptide (CNP) (1-22), human
Cat No. : RP11110
Lot No. : E9189224K
Time Processed :13:33:36
Month-Day-Year Processed :11/19/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm
Equipment: ZJ19010324

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	5.617	68904	9991	1.294
2	10.883	14695	2111	0.276
3	11.175	50950	7063	0.957
4	11.383	5059155	555894	95.000
5	12.658	17592	2628	0.330
6	20.175	114103	12063	2.143
Total		5325399	589750	100.000