


Mass Spectrum



Sample Information

Month-Day Processed : 11/13/24
Time Processed : 19:51:18
Injection Volume : 0.2

Sample Name : C-Type Natriuretic Peptide (CNP) (1-22), human

Cat No. : RP11110

Lot No. : E9189224K

Theoretical MW: 2197.60

Observed MW : 2197.4

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : ZJ21010035
Interface Bias : -3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH