

Sample Name: Melittin,honey bee
Cat No. : RP20415
Lot No. : E9207424I
Time Processed: 16:04:58
Month-Day-Year Processed :09/28/2024
Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm
<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	Pump A B.Conc	5
25.00	Pumps	Pump A B.Conc	65
25.01	Pumps	Pump A B.Conc	95
27.00	Pumps	Pump A B.Conc	95
27.01	Pumps	Pump A B.Conc	5
35.00	Pumps	Pump A B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column : Inertsil ODS-SP 4.6 x 250 mm

Equipment:ZJ19010014

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	17.400	11024	1417	0.072
2	17.748	59983	2956	0.392
3	18.625	52992	1626	0.347
4	18.750	7753	1558	0.051
5	18.913	14752	1565	0.097
6	19.333	51353	2092	0.336
7	19.511	20117	2020	0.132
8	20.043	263848	14355	1.726
9	20.692	197012	13164	1.289
10	21.153	14525390	944138	95.034
11	22.208	5877	1200	0.038
12	22.300	5604	1129	0.037
13	22.389	6863	1091	0.045
14	22.649	18245	1131	0.119
15	22.775	4431	998	0.029
16	23.011	26225	1027	0.172

Peak#	Ret. Time	Area	Height	Area%
17	23.469	5089	408	0.033
18	23.635	2117	304	0.014
19	23.712	1447	282	0.009
20	23.820	1285	224	0.008
21	23.953	1576	209	0.010
22	24.101	1501	176	0.010
Total		15284483	993067	100.000