


Mass Spectrum



Sample Information

Month-Day Processed :	09/29/22	Equipment	: ZJ21010035
Time Processed :	16:31:42	Interface Bias	: +4.5 kV
Injection Volume :	0.4	Drying Gas Flow	: 5 L/min
Sample Name :	DYKDDDDK Peptide	T.Flow	: 0.2 ml/min
Cat No. :	RP10586-1	B.conc	: 50% H2O/50% MeOH
Lot No. :	E9259422I		
Theoretical MW :	1012.97		
Observed MW :	1012.4		