

Sample Name :OVA Peptide (323-339)
Cat No. :RP10610
Lot No. :E9292126H

Time Processed :14:50:35

Month-Day-Year Processed :08/23/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ19010324

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	5.125	11693	1933	0.302
2	6.925	3684544	543096	95.003
3	7.108	121019	23138	3.120
4	7.375	20781	2979	0.536
5	7.550	30748	2904	0.793
6	7.717	9566	1552	0.247
Total		3878352	575601	100.000