


Mass Spectrum



Sample Information

Month-Day Processed :	03/24/23	Equipment	: ZJ21010035
Time Processed :	7:51:26	Interface Bias	: +4.5 kV
Injection Volume :	0.3	Drying Gas Flow	: 5 L/min
Sample Name :	DYKDDDDDK Peptide	T.Flow	: 0.2 ml/min
Cat No. :	RP10586	B.conc	: 50% H2O/50% MeOH
Lot No. :	E9337623C		
Theoretical MW: 1012.97			
Observed MW : 1012.8			