

Sample Name: Influenza A NP(366-374) Strain A/PR/8/35

Cat No. : RP20267

Lot No. : E9400724C

Time Processed : 17:40:09

Month-Day-Year Processed : 03/20/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow: 1 ml/min

Wavelength: 220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column : Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ19010140

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	6.608	172953	31752	2.434
2	7.017	30415	5620	0.428
3	9.067	6819425	821718	95.981
4	9.517	21801	3126	0.307
5	9.592	9990	2002	0.141
6	10.092	38291	6060	0.539
7	10.242	12110	2149	0.170
Total		7104985	872428	100.000