

Sample Name: Influenza A NP(366-374) Strain A/PR/8/35

Cat No. : RP20267

Lot No. : E9433825B

Time Processed : 8:51:25

Month-Day-Year Processed : 02/08/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow: 1 ml/min

Wavelength: 220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm

Equipment: ZJ24010222

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	7.775	8635	1365	0.362
2	8.250	2335997	374229	98.052
3	12.667	11112	1009	0.466
4	13.533	14942	2149	0.627
5	14.558	11717	1763	0.492
Total		2382403	380515	100.000