

Sample Name :Influenza A NP(366-374) Strain A/PR/8/35

Cat No. : RP20267

Lot No. : E9523724H

Time Processed :6:14:35

Month-Day-Year Processed :08/10/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: GK12010012

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	7.183	5509857	652771	98.403
2	7.550	26708	3832	0.477
3	7.708	17985	1931	0.321
4	8.283	22883	2766	0.409
5	8.483	21862	1963	0.390
Total		5599295	663262	100.000