

Sample Name :TAT peptide
 Cat No. :RP20256
 Lot No. :E9607526E
 Time Processed :11:28:17
 Month-Day-Year Processed :05/27/2026
 Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm
 <<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	0
5.00	Pumps	B.Conc	0
20.00	Pumps	B.Conc	30
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	0
35.00	Pumps	B.Conc	0
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ25010054

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	11.217	11990699	1814489	97.001
2	11.583	142196	12436	1.150
3	11.900	135109	13931	1.093
4	12.392	36396	3074	0.294
5	12.567	12377	1518	0.100
6	24.467	44699	1441	0.362
Total		12361475	1846889	100.000