


# Mass Spectrum



## Sample Information

|                          |                                                |                 |                 |
|--------------------------|------------------------------------------------|-----------------|-----------------|
| Month-Day Processed :    | 07/26/26                                       | Equipment       | :Z126010077     |
| Time Processed :         | 8:52:25                                        | Interface Bias  | : +4.5 kV       |
| Injection Volume :       | 0.5                                            | Drying Gas Flow | :5 L/min        |
| Sample Name :            | C-Type Natriuretic Peptide (CNP) (1-22), human | T.Flow          | :0.2 ml/min     |
| Cat No. :RP11110         |                                                | B.conc          | :50%H2O/50%MeOH |
| Lot No. :E9659426G       |                                                |                 |                 |
| Theoretical MW : 2197.60 |                                                |                 |                 |
| Observed MW : 2197.8     |                                                |                 |                 |