

Sample Name :Tau peptide (275-305) (Repeat2 domain)

Cat No. : RP30887

Lot No. : E9665125E

Time Processed :4:06:48

Month-Day-Year Processed :05/19/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ20010140

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	10.742	9092644	966111	95.208
2	10.933	326794	82002	3.422
3	11.117	26900	7164	0.282
4	12.417	103992	16292	1.089
Total		9550329	1071568	100.000