


# Mass Spectrum



## Sample Information

Month-Day Processed : 05/19/25  
 Time Processed : 8:07:49  
 Injection Volume : 0.2  
 Sample Name : Tau peptide (45-73) (Exon2/Insert1 domain)  
 Cat No. : RP30888  
 Lot No. : E9666025E  
 Theoretical MW : 2977.99  
 Observed MW : 2977.2

## Interface

Interface : ESI  
 Nebulizing Gas Flow : 1.5L/min  
 CDL Temp : 250  
 Block Temp : 200

## Equipment

Equipment : ZJ21010035  
 Interface Bias : -3.5 kV  
 Drying Gas Flow : 5 L/min  
 T.Flow : 0.2 ml/min  
 B.conc : 50%H2O/50%MeOH