

Sample Name :Melittin,honey bee
Cat No. :RP20415
Lot No. : E9674425H
Time Processed :2:47:52
Month-Day-Year Processed :08/28/2025
Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ17010504

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	10.233	56913	7663	0.997
2	18.117	23335	3204	0.409
3	18.475	5422515	507952	95.018
4	18.867	145584	11759	2.551
5	19.050	37270	4257	0.653
6	19.292	21228	1290	0.372
Total		5706844	536124	100.000