


Mass Spectrum



Sample Information

Month-Day Processed : 05/18/25
 Time Processed : 22:13:34
 Injection Volume : 0.2
 Sample Name : Neurogranin (48-76), mouse
 Cat No. : RP30891
 Lot No. : E9677825E
 Theoretical MW: 2322.46
 Observed MW : 2322.4

Interface :ESI
 Nebulizing Gas Flow :1.5L/min
 CDL Temp :250
 Block Temp :200

Equipment
 Interface Bias
 Drying Gas Flow
 T.Flow
 B.conc

:ZJ22010150
 :+4.5 kV
 :5 L/min
 :0.2 ml/min
 :50%H₂O/50%MeOH