


Mass Spectrum



Sample Information

Month-Day Processed :	07/07/22	Equipment	: ZJ21010035
Time Processed :	19:45:32	Interface Bias	: +4.5 kV
Injection Volume :	0.4	Drying Gas Flow	:5 L/min
Sample Name :	TAT peptide	T.Flow	:0.2 ml/min
Cat No. : RP20256		B.conc	:50%H2O/50%MeOH
Lot No. : E9685922F			
Theoretical MW: 1621.91			
Observed MW : 1621.6			