

Sample Name :[Gly22]- β -Amyloid (1-42), Arctic Mutation

Cat No. : RP30247

Lot No. : E9704025H

Time Processed :17:32:06

Month-Day-Year Processed :09/11/2025

Pump A : 0.1% Ammonia solution in 100% water (v/v)

Pump B : 0.1% Ammonia solution in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :UniPS 10-300 10um,4.6x250mm

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	9.600	117364	8242	0.944
2	10.267	342742	20453	2.757
3	10.992	11823589	459428	95.096
4	12.108	149576	5927	1.203
Total		12433271	494050	100.000