


Mass Spectrum



Sample Information

Month-Day Processed :	04/26/25	Interface :	ESI	Equipment	:ZJ22010150
Time Processed :	16:20:19	Nebulizing Gas Flow :	1.5L/min	Interface Bias	: +4.5 kV
Injection Volume :	0.2	CDL Temp :	250	Drying Gas Flow	:5 L/min
Sample Name :	C-Type Natriuretic Peptide (CNP) (1-22), human	Block Temp :	200	T.Flow	:0.2 ml/min
Cat No. :	RP11110			B.conc	:50%H2O/50%MeOH
Lot No. :	E9705525D				
Theoretical MW :	2197.60				
Observed MW :	2197.5				